

SIEMENS



ADVIA Centaur® XPT System
PEP to Go

Answers for life.

Table of Contents

1. Hardware Overview
2. Software Overview
3. Loading Supplies and Reagents
4. Sample Processing
5. Monitoring Samples and Results
6. Calibration
7. Quality Control
8. Daily Maintenance
9. Weekly Maintenance
10. Monthly Maintenance
11. As Needed Maintenance
12. Troubleshooting
13. System Restart



Hardware Overview

- A. Sample Tip Loader
- B. Cuvette Loader
- C. Sample Entry Queue
- D. STAT Entry
- E. Sample Exit Queue
- F. Ancillary Reagent Entry
 - 25 positions in Ancillary Compartment
- G. Cleaning Solution and Sample Tip Waste
- H. System Fluids
 - Acid, Base, and Wash 1
- I. Tip Tray Waste and Cuvette Waste
- J. Waste/Water Drawer
- K. Primary Reagent Compartment
 - 30 positions



Software Overview

A. Command Bar


Provides access to all system functions and information required to operate the system

B. Workspace

Displays a set of tabs or a window containing tasks or information when a Command Bar button is selected

C. Status bar


Provides access to windows that display information about the entire system



Software Overview

Status Bar

- A. Online Help
- B. Sign-In with Operator Name and Password
- C. Change the System State
- D. File Explorer
- E. Remote Connection
- F. Change the status of the LIS connection
- G. LAS connection – view status and reset
- H. Dashboard – quick count of all samples and tests currently managed by the system
- I. STAT Samples screen
- J. Find Patient/QC results
- K. System Status Display
- L. Turn on/off alarm sounds



Loading Supplies and Reagents

Loading Supplies and Emptying Waste

Cuvettes

- May be loaded at any time
- The cuvette loading bin holds up to 1000 cuvettes (approximately 5 bags)

Sample Tips

- May be loaded at any time
- Each tip tray holds 120 tips. Up to 840 tips can be loaded onto the system

Deionized Water Bottle

- May be removed from the system at any time
- The capacity of the water bottle is 10L
- Option for direct plumbing

Acid, Base, and Wash 1 Bottles

- May be replaced at any time
- Should be replaced when any of the following conditions exist:
 - They have been on the system for a month
 - The reagent has expired
 - The reagent has insufficient volume

Waste

- The cuvette waste, tip waste, tip tray waste, and liquid waste may be emptied at any time
- After emptying waste, respond appropriately to the popup on the Workspace
- Option for liquid waste direct plumbing



Note:

Replace acid and base bottle at the same time to manage lot numbers and to obtain optimum reagent performance. It is important to place the acid and base bottles in the correct positions on the system.

Loading Primary Reagents

1. Mix by hand to resuspend the paramagnetic particles of the Solid Phase reagent
 - Raise one end of the pack to a 90° angle and return to a horizontal position
 - Raise the other end of the pack to a 90° angle and return to a horizontal position
 - Repeat at least 20 times or until the paramagnetic particles are no longer visible on the bottom of the pack
2. Load the pack
 - Open the primary reagent compartment and identify a position for the pack
 - Pull the reagent holder toward you and load the pack
 - Push the reagent holder back
 - Close the primary reagent compartment door


Note:

For optimal throughput, install packs of similar types next to one another. Similar types of reagents are identified by the color on the end of the label. Certain methods require special primary reagent handling. Refer to Instructions for Use for more details.



Removing Primary Reagents

1. At the Command Bar, select **Reagents**
2. The first tab is **Primary Reagent Compartment**
3. Select the reagent pack
4. Select **Remove**
5. After the system verifies the pack is no longer in use, the Safe to Unload indicator displays next to the pack



Note:

Do not remove a primary pack when the system is in process if the LED on the primary reagent door is green.


Loading Ancillary Reagents

1. Gently invert ancillary reagent packs several times while minimizing foaming
2. Place the reagent pack in the ancillary reagent entry with the barcoded side facing outward
3. An LED under the ancillary reagent entry turns green to indicate that the pack is loading
4. The system moves the pack into the ancillary reagent compartment



Removing Ancillary Reagents


1. At the Command Bar, select **Reagents**
2. Select the **Ancillary Reagent Compartment** tab
3. Select the ancillary pack
4. Select **Remove**
5. The system ejects the pack into the ancillary reagent entry
6. An LED next to the ancillary reagent entry turns orange to indicate that the pack is exiting the system



Sample Processing

Manually Ordering Patients By SID


1. At the command bar, select **Orders**
2. The first tab is **Create Patient Orders**
3. Enter the Sample Identification Number (SID) and then select the Enter key on the keyboard to activate the test menu
4. In the Sample Information area, select a Specimen Type and other sample information, as needed
5. In the Tests area, select one or more tests. Test profiles and dilution profiles can also be selected
6. Select the **Save** button



Entering SID for Non-Barcoded Samples

For samples onboard the system with test orders:


1. At the command bar, select **Samples**
2. The first tab is **InProcess Queue**
3. Select the rack that contains the unidentified sample tube
4. On the Rack and Tube diagram, select the unidentified sample tube
5. Select **Enter Barcode**
6. Select the type of sample
7. Enter the barcode for the sample in the Sample ID text box
8. Select the **OK** button



Processing STATs


There are three options for processing STAT patient samples:

1. LIS – STAT orders received from the LIS are given priority over routine samples
2. STAT Entry – load racks in the STAT entry that require processing before racks in the sample entry queue
3. Create Patient Orders Screen – select the STAT priority when ordering the sample




Ordering Dilutions and Replicates

1. At the command bar, select **Orders**
2. The first tab is **Create Patient Orders**
3. Enter the SID and then select the Enter key on the keyboard
4. In the Sample Information area, select a Specimen Type and other sample information, as needed
5. For a Manual Dilution (prepared on the bench): In the Manual Dilution field, enter the manual dilution factor
6. In the Tests area, select the test
7. For a System Performed Dilution (performed by the system using an onboard diluent): Under the Dilution heading, select the appropriate dilution factor
8. If needed, select the number of replicates and reagent pack information. If there are multiple reagent lots onboard, the reagent lot can be selected
9. Select the **Save** button



Ordering Batches


1. At the command bar, select **Orders**
2. Select the **Create Batch Orders** tab
3. Enter a starting rack number and an ending rack number
4. Select a specimen type
5. Enter a manual dilution factor and comments if applicable
6. Select the tests to run. The system performs the selected tests on all samples loaded in the starting rack, ending rack and all racks in between
7. Select the **Save** button
8. To view the status of the batch order, select the **View Batch Orders** tab



Monitoring Samples and Results

View Pending Orders screen

- To access, select **Orders** from the Command Bar, and then select the **View Pending Orders** tab
- Used to view QC and patient test orders that have not yet started processing
- Options are available to:
 1. Edit orders
 2. Delete orders
 3. Query LIS



Test Results Overview screen

- To access, select **Test Results** from the Command Bar, and then select the **Overview** tab
- Used to view test orders and results
- Options are available to:
 - Post an order to the Watch List screen – used to selectively monitor certain test orders
 - Repeat a test
 - Transmit results to the LIS
 - Export data
 - Print reports

The screenshot displays the Siemens test results overview screen. The interface includes a top navigation bar with icons for Samples, Reports, Calibration, Supplies, Maintenance, Diagnostics, Inprocess, Orders, Test Results, QC, Setup, Utilities, and Events. The 'Test Results' tab is active, and the 'Overview' sub-tab is selected. A search bar at the top right shows 'Find in Page: SD 12345' and 'Page 1 of 1. Total tests: 23'. Below the search bar, there are filters for 'Today's Samples' and 'Orders Created or Updated Today'. The main area contains a table of test orders with columns for Order ID, Date & Time, Test, Result, and Interpretation. The table lists several orders with their respective results and interpretations. On the right side, there are buttons for 'Sample', 'Edit/View...', 'Start Watching', 'Delete', 'Move to Historical', 'Test', 'Accept', 'Repeat', 'Assign Cal', 'Transmit', and 'Cancel'. At the bottom, there are buttons for 'Export', 'Print...', and 'Close'. Five orange circles with numbers 1 through 5 are overlaid on the screen to indicate key features: 1 points to the 'Start Watching' button, 2 points to the 'Repeat' button, 3 points to the 'Transmit' button, 4 points to the 'Export' button, and 5 points to the 'Print...' button.

Order ID	Date & Time	Test	Result	Interpretation
1220	3/7/2014 6:33:56 AM	ES	+ 2.53 ng/mL	Above Linearity, > Dilution ...
		TA	RU ng/dL	Cal Invalid
		TA	+ 8.00 ng/mL	> Conc Range, Autorepeat
		TA	+ 12.7 ug/dL	> Conc Range, Autorepeat
1127	3/7/2014 6:34:27 AM	TA	+ 8.00 ng/mL	> Conc Range, Autorepeat
		TA	+ 13.1 ug/dL	> Conc Range, Autorepeat
1126	3/7/2014 6:34:21 AM	TA	+ 8.00 ng/mL	> Conc Range, Autorepeat
		TA	+ 30.0 ug/dL	> Conc Range
1125	3/7/2014 6:34:16 AM	TA	+ 8.00 ng/mL	> Conc Range, Autorepeat
		TA	+ 30.0 ug/dL	> Conc Range

Dashboard screen

- To access, select the Dashboard icon from the Status Bar at the bottom of the screen
- Used to post a quick view of samples and tests by their current state
- The system updates the counters for the samples and the tests in real time until the sample or test is deleted or moved to Historical
- Selecting one or more counter on either tab launches the Test Results Dashboard tab. The counters selected will be highlighted

The screenshot displays the Siemens software interface for monitoring sample and test status. The top navigation bar includes icons for Samples, Requests, Calibration, Supplies, Maintenance, Diagnostics, Inprocess, Orders, Test Results, QC, Setup, Utilities, and Events. The 'Test Results' tab is active, showing a 'Dashboard: Patient and QC' window. This window features a 'Sample Status: Inprocess' table with columns for 'Order Date & Time', 'Test', and 'Result'. A modal window is open over the table, showing a summary of test counts: Pending test (7), Inprocess (1), Intervention Needed (0), Completed (18), and Transmitted (0). The bottom status bar contains various system icons, with the 'Dashboard' icon highlighted by a yellow box. The system date and time are displayed as 'Friday, March 07, 2014 8:25:32 AM'.

Order ID	Order Date & Time	Test	Result
22225	3/7/2014 8:23:26 AM	OS	2
22225	3/7/2014 8:24:07 AM	TS	2
22224	3/7/2014 8:23:16 AM	OS	2
22224	3/7/2014 8:23:16 AM	TS	2
22223	3/7/2014 8:24:06 AM	OS	2
22223	3/7/2014 8:23:08 AM	TS	2
12347	3/7/2014 8:07:33 AM		

Calibration

Reviewing Calibration Status

1. At the Command Bar, select **Reagents**
2. Select the **Status** tab
3. Review the calibration status and expiration date for the onboard reagents

The screenshot shows the Siemens reagent status interface. The Command Bar at the top has the 'Reagents' tab selected. The 'Status' tab is also selected. The table below shows the following data:

Reagent Type	Reagent Name	Lit.	Position	Pack Status	Test Name	Test Status	Inventory	Lot Expiration	Remaining DBI	Calibration Status	Calibration Expiration	Barcode
Primary Compartment												
	Primary CoA	017	12		CoA	100 Tests	4/12/2015	424 Db		Not Calibrated		100171000001
	Primary DIG	218	11		DIG	100 Tests	10/9/2014	424 Db		Current Calibration	3/19/2014 10:18:27 AM	0212101000001
	Primary FT4	047	1		FT4	100 Tests	9/20/2014	424 Db		Expired	2/27/2014 1:25:18 PM	1136471000001
	Primary RuBtA	157	21		RuBtA	100 Tests	4/11/2014	284 Db		Not Calibrated		0841571000001
	Primary T3	184	3		T3	100 Tests	8/15/2014	424 Db		Current Calibration	3/27/2014 8:25:40 AM	0041841000001
	Primary T4	109	3		T4	100 Tests	5/6/2014	424 Db		Current Calibration	3/11/2014 11:28:54 AM	0011910000001
	Primary TnCG	204	4		TnCG	100 Tests	8/9/2014	216 Db		Expired	3/28/2014 7:55:18 AM	0292410000001
	Primary TSH	261	5		TSH	100 Tests	7/21/2014	216 Db		Current Calibration	3/26/2014 8:29:07 AM	0012811000001
Auxiliary Compartment												
	Diluent ADI4.2	254	3			10,000 mL	8/13/2014					0911254000001
	Diluent ADI4.4	1674	4			5,000 mL	6/18/2014					0940167400001
	Auxiliary Wash T3F4B12	2804	1			25,000 mL	10/7/2014					0030280400001
	Auxiliary T3TOR	3214	2			5,000 mL	11/17/2014					0363214000001

Note: Calibration status can also be viewed on the Calibration Summary screen. Select **Calibration** from the Command Bar, and then select the **Calibration Summary** tab.

Scanning Master Curve Test Definition (MC DEF) Cards

1. Locate the MC DEF card in the primary reagent box. Select the ADVIA Centaur® side of the card
2. Scan the first 2D barcode on the MC DEF card by placing the card on the barcode reader platform to the right of the touchscreen. Pass the card under the scanner until the blue light points at the center of the 2D barcode. The scanner beeps once and the green light on the scanner flashes once after successfully reading the first 2D barcode
3. Scan the second 2D barcode. The scanner beeps twice and the green light on the scanner flashes twice after successfully reading the second 2D barcode
4. Scan the third and fourth 2D barcodes if included on the MCDEF card. The scanner beeps 3 times and 4 times after successfully reading the third and fourth 2D barcodes
5. Select **Save** in the confirmation window to save the master curve and test definition and close the window



Note:

The MC DEF cards must be scanned with every new lot of reagent. The MC DEF cards can only be scanned when the system is in the Ready state.


Scanning Calibrator Definition Cards

1. Locate the Calibrator Definition card in the calibrator kit. It is labeled ADVIA Centaur® Assigned Value Card and has a purple bar at the top
2. At the Command Bar, select **Calibration**
3. Select the **Calibrator Definitions** tab
4. Select the **Add Calibrator Definition** button
5. Scan all of the barcodes on the Calibrator Definition card from top to bottom
6. Select the **Save** button

Note:

The Calibrator Definition cards must be scanned with every new lot of calibrator.

The Calibrator Definition cards can be scanned while the system is In Process.



Manually Ordering Calibrations

1. At the Command Bar, select **Calibration**
2. Select the **Calibration Order** tab
3. Select a test to calibrate
4. Select a reagent lot
5. Select a calibrator lot
6. Select the **Save** button

SIEMENS IDEMA Control

Ready

Calibration Order

Test

Reagent Lot

Select a calibrator lot for DIG, reagent lot 209

Lot	Expiration Date
209	6/18/2014
210	10/9/2014

Name	Lot	Expiration Date
CALB	CB40	4/9/2015

Summary

Test	Reagent Lot	Calibrator Name	Calibrator Lot	Priority	Replicates	Reagents	Level	Sample ID	Back	Position
DIG	209	CALB	CB40	Routine	1	P...	Low	CB40		
							High	CB40H		

Save Clear Done

Friday, March 07, 2014 8:12:16 AM

Reviewing Calibration Data

1. At the Command Bar, select **Calibration**
2. Select the **Calibration Results Overview** tab
3. Review the Result column for each calibration test
 - A Valid calibration indicates that the Calibration Evaluation Criteria have passed
 - An Invalid calibration indicates that one of the criteria did not pass and further investigation will be required
4. To view more information on an Invalid calibration, select the test and then select the **Calibration Data** button

The screenshot displays the Siemens ADMVA Control System interface. The Command Bar at the top shows the 'Calibration' icon selected. Below it, the 'Calibration Results Overview' tab is active. The main area contains a table of calibration tests with columns for Test, Reagent Lot, Calibrator Lot, Calibration Ordered, Order Status, Result Completed, and Result. A 'Calibration Data' button is located on the right side of the table. The interface also includes a 'Ready' status indicator, various system icons, and a footer with system state and date information.

Test	Reagent Lot	Calibrator Lot	Calibration Ordered	Order Status	Result Completed	Result
Coef	021	C201	3/4/2014	Historical		Valid
DIG	210	C040	2/19/2014	Historical		Invalid
DIG	210	C040	2/19/2014	Historical		Valid
EHV	152	C152	2/27/2014	Historical		Invalid
EHV	152	C152	2/27/2014	Historical		Valid
FiskeIDA	225	C725	2/19/2014	Historical		Valid
FiskeIDA	225	C725	2/27/2014	Historical		Valid
FTx	047	C475	2/28/2014	Historical		Valid
HA	011	C104	2/28/2014	Historical		Valid
HA/IT	055	C145	2/21/2014	Historical	Inst Integrity, System Error	
HA/IT	055	C145	2/21/2014	Historical		Valid
HbCT	061	C501	2/14/2014	Historical		Invalid
HbCT	061	C501	2/17/2014	Historical		Invalid
HbCT	061	C501	2/17/2014	Historical		Invalid
HbCT	061	C501	2/17/2014	Historical		Invalid
HbCT	061	C501	2/17/2014	Historical		Invalid

Accepting a Calibration

- Required when the calibration status is "Pending Accept" status on the Calibration Results Overview screen
- "Pending Accept" calibrations are outside of observed ranges but within defined ranges
 - Observed ranges: After 4 calibrations of the same reagent lot and calibrator lot combination, the system calculates observed ranges
 - Defined ranges: The values for Slope, Ratio, Low Cal Deviation and High Cal Deviation entered from the MC DEF card
- To accept, select the **Accept** button on the Calibration Data screen
- After accepting, the calibration status changes to "Valid Operator Accepted"

The screenshot displays the Siemens ADMA Control software interface. The top navigation bar includes icons for Samples, Reagents, Calibration, Supplies, Maintenance, Diagnostics, and a 'Ready' status indicator. The main window is titled 'Calibration Data' and shows the following information:

- Calibration Status:** Current Calibration, Test: TSH, Calibrator Lot: CB4B, Minimum Replicates: 2
- Predefined Filters:** Current: Editable, Reagent Lot: 261, Calibrator Name: CALB, Critical Precision: 0.255
- Calibration History Table:**

Calibration	Calibration Result	Related Calibration	Calibrator Lot
3/7/2014 5:52:18 AM	Pending Accept		CB4B
3/7/2014 9:48:27 AM	Valid		CB4B
- Calibrated:** 3/7/2014 5:52:18 AM
- Evaluation Ranges Table:**

Range Type	Calibration Value	Defined Ranges	Observed Ranges
Slope	0.826	0.330 - 4.670	0.685 - 1.185
Ratio	407.586	134.719 - 1283.290	425.700 - 771.726
Low Cal Deviation	0.357	-1.070 - 0.670	0.060 - 0.263
High Cal Deviation	0.135	-1.070 - 0.670	-0.880 - 0.176
- Calibrator Target Values and Observed Results Table:**

Level	ID	Conc (µg/mL)	Expected RLU	Replicate Number	Measured RLU	Acceptable CV %	Calculated CV %
High	CB4BH	32.90	84028	1	102454.50	8.2	0.128
				2	1023456		
Low	CB4BL	0.00	1171	1	2182.00	19.2	11.018
				1	2352		
- Buttons:** Master Curve Definitions, View, Plot, Close

The 'Accept' button on the right side of the 'Calibration Data' window is highlighted with a yellow box.

Quality Control

Adding a QC Definition

1. At the Command Bar, select **QC**
2. Select the **QC Definitions** tab
3. Select the **Add QC Definition** button. The Add New QC Definition screen displays
4. Enter the appropriate information in the Control Details area
5. Select a test
6. Enter the expected values – mean and 2SD or low and high range
7. Add additional tests as needed
8. Select the **Save** button

Editing a QC Definition

To edit an existing QC definition, on the QC Definition screen, select the QC definition to edit and then select the **Edit/View** button. Make the necessary edits and select the **Save** button.

The top screenshot shows the Siemens software interface with the 'QC Definitions' tab selected. A table lists various QC definitions with columns for Control Name, Control Level, Control ID, Control Lot, Control Lot Expiration, Stage, and Date Last Used. A red circle '1' is placed over the 'QC' icon in the command bar, and a red circle '2' is placed over the 'QC Definitions' tab.

The bottom screenshot shows the 'Add New QC Definition' screen. It has several fields: Control Name (QC1), Control ID (QC1), Control Lot (12345), Control Lot Expiration (4/30/2015), Stage (Primary), and Control Type (Routine). A red circle '4' is placed over the Control Name field. Below these fields is a list of tests: XPO, CA1B, DH, COB, CA, DG (checked), ELP, and FusedA. A red circle '5' is placed over the DG test. To the right is a table for Test Name, Units, Low, High, and Target (Mean) Range. A red circle '6' is placed over the Target Range field. At the bottom right are 'Save' and 'Cancel' buttons, with a red circle '8' placed over the 'Save' button.

Manually Ordering QC

1. At the Command Bar, select **Orders**
2. Select the **QC Orders** tab
3. Select a priority
4. In the Tests area, select one or more tests
5. To order a QC profile, select the **Select Profiles** button and then select the desired profile
6. In the Control Materials area, select each control material to include in this order
7. Select the **Add** button
8. In the Orders area, select or enter the order and test information as needed by selecting the Details/Summary button
9. Select the **Save** button

The screenshot shows the Siemens ADVA Control software interface. The top command bar includes buttons for Samples, Reagents, Calibration, Supplies, Maintenance, Diagnostics, Ready, Orders, Test Results, QC, Setup, Utilities, and Events. The 'Orders' tab is selected. Below the command bar, there are tabs for 'Create Patient Orders', 'Create QC Orders', 'View Pending Orders', 'Create Batch Orders', 'View Batch Orders', and 'Control Buckets'. The 'Create QC Orders' tab is active. The 'Priority' section has 'Routine' selected. The 'Tests' area shows a list of tests: DIG, FTA, GENI, T1, T4, ThCG (checked), and VBIG. The 'Control Materials' area is empty. The 'Orders' table at the bottom contains three rows of order information:

Control Name	Level	Control ID	Control Lot	Lot Expiration	Usage	Rack ID	Test Name	RegCofact	Reagents
CONTROLPLLS1...	1	1001	4796	2/26/2015	Primary		ThCG	1	Lot#: Pack ID: Ancillary: Ancillary2: Diluent: onlyReagentLotProvided:False
CONTROLPLLS1...	2	1002	1232	2/19/2015	Primary		ThCG	1	Lot#: Pack ID: Ancillary: Ancillary2: Diluent: onlyReagentLotProvided:False
CONTROLPLLS1...	3	1003	1111	2/26/2015	Primary		ThCG	1	Lot#: Pack ID: Ancillary: Ancillary2: Diluent: onlyReagentLotProvided:False

The bottom of the interface has a 'System State' button and a status bar showing 'Friday, March 07, 2014 1:30:53 AM'.

Printing QC Data

1. At the Command Bar, select **QC**
2. Select the **QC Statistics** tab
3. Select either the **Review** or **Analysis** tab
4. Select **Reports**

The screenshot displays the SIEMENS software interface for QC Statistics. The interface includes a Command Bar at the top with icons for Reagents, Calibration, Supplies, Maintenance, and Diagnostics. The QC icon is highlighted with a red circle 1. Below the Command Bar, the QC Statistics tab is selected, highlighted with a red circle 2. The Review and Analysis tabs are visible, with the Review tab highlighted by a red circle 3. The main area shows a table of test results with columns for Date, Analyzer, Module, Test, Control, Level, Lot, QC Violation, Result, Unit, Z, Target ± 2SD, Comment, Excluded, Reviewed, and Reviewed by. The table contains 15 rows of data, with several rows highlighted in red. A red circle 4 highlights the Reports button at the bottom right of the table. The bottom status bar shows the system state and date/time: Friday, March 07, 2014 12:42:34 AM.

Date	Analyzer	Module	Test	Control	Level	Lot	QC Violation	Result	Unit	Z	Target ± 2SD	Comment	Excluded	Reviewed	Reviewed by
3/18/2014 9:44:25...	EP-102	FTx	Control 2	2	65422			2.28	ng/dL		-1.7 2 - 4.5			X	System
3/18/2014 9:44:19...	EP-102	FTx	Control 1	1	54221			0.76	ng/dL		-1.5 0.5 - 2.5			X	System
3/17/2014 9:53:06...	EP-102	FTx	Control 3	3	76543			6.37	ng/dL		1.2 4 - 7			X	System
3/17/2014 9:32:58...	EP-102	FTx	Control 2	2	65422			2.29	ng/dL		-1.7 2 - 4.5			X	System
3/17/2014 9:32:54...	EP-102	FTx	Control 1	1	54221			0.76	ng/dL		-1.5 0.5 - 2.5			X	System
3/13/2014 3:07:18...	EP-102	FTx	Control 3	3	76543			6.37	ng/dL		1.2 4 - 7			X	System
3/13/2014 3:06:54...	EP-102	FTx	Control 2	2	65422		[N] [3]	5.53	ng/dL		3.4 2 - 4.5		X	X	labmanag
3/13/2014 3:06:38...	EP-102	FTx	Control 1	1	54221			0.76	ng/dL		-1.5 0.5 - 2.5			X	System
3/13/2014 10:42:11...	EP-102	FTx	Control 2	2	65422		[N] [3]	5.42	ng/dL		4.2 2 - 4.5		X	X	labmanag
3/12/2014 10:09:4...	EP-102	FTx	Control 3	3	76543			6.37	ng/dL		1.2 4 - 7			X	System
3/12/2014 10:09:0...	EP-102	FTx	Control 2	2	65422		[N] [3]	4.88	ng/dL		5.8 2 - 4.5		X	X	labmanag
3/12/2014 8:56:24...	EP-102	FTx	Control 2	2	76543		[N] [3]	2.28	ng/dL		-4.3 4 - 7			X	labmanag
3/12/2014 8:56:08...	EP-102	FTx	Control 2	2	65422		[N] [3]	6.41	ng/dL		5.4 2 - 4.5		X	X	labmanag
3/12/2014 8:55:43...	EP-102	FTx	Control 1	1	54221			0.75	ng/dL		-1.5 0.5 - 2.5			X	System
3/10/2014 3:59:38...	EP-102	FTx	Control 1	1	54221			0.71	ng/dL		-1.6 0.5 - 2.5			X	System

Exporting QC Data

1. At the Command Bar, select **QC**
2. Select the **QC Statistics** tab
3. Select the **Tools** tab
4. Select the **Export** tab
5. Select the data to export, the time filter, and the file type
6. Select **Export**


The screenshot shows the Siemens software interface with the following elements:

- Command Bar:** Includes icons for Samples, Reagents, Calibration, Supplies, Maintenance, Diagnostics, Ready, Orders, Test Results, QC (highlighted with a red circle 1), Setup, Utilities, and Events.
- QC Statistics Tab:** Selected (highlighted with a red circle 2). It contains sub-tabs: Review (highlighted with a red circle 4), Messages, Tests, Setup, Maintenance, Export (highlighted with a red circle 3), and Database Backup.
- Table:** A table with columns: Test, Severity, Level, Lot, File Size, Analyzer, Module, Unit. It lists tests: AFP, DGC (selected), and FT4.
- Export Dialog:** A dialog box titled 'Export' with a 'Show by' dropdown set to 'Test'. It lists items to export: Test (selected), Control, Lot, Result, Date, Excluded, Range Flag, System Flag, User, and Z. The 'Export' button is highlighted with a red circle 6.
- Status Bar:** Shows 'Friday, March 07, 2014 11:03 AM'.

Daily Maintenance

Daily Cleaning Procedure

1. Check/replenish cleaning solution (7 day onboard stability)
 - Add 1 container of ADVIA Centaur® Cleaning Solution Concentrate to the cleaning solution bottle
 - Fill the cleaning solution bottle to the neck with DI water
 - Place the bottle onto the system
2. Check/replenish cuvettes
3. Check/replenish the water bottle
4. Check/empty the liquid waste bottle
5. Check/empty the cuvette waste bin
6. Start the Daily Cleaning Procedure:
 - A. Select **Maintenance** on the Command Bar
 - B. Select **Daily Maintenance** on the Schedule tab
 - C. Select **Perform**



Weekly Maintenance

Cleaning the Water Bottle

1. Remove the cap and empty any water from the water bottle
2. Pour approximately 1L of cleaning solution into the water bottle
3. Place the cap on the bottle and close tightly
4. Swirl the solution around the inside of the bottle a few times
5. Press some gauze over the holes in the top of the bottle
6. Invert the bottle, swirl the solution a few times and then return the bottle to the upright position
7. Soak the bottle in the upright position for 5 minutes
8. Repeat steps 5 and 6
9. Swirl the solution again and then empty the bottle
10. Rinse the inside of the bottle including the tubing at least 5 times with fresh reagent water
11. Invert the bottle to dry



Note:

This procedure is not required if the system is direct plumbed for water.

Emptying the Water Trap

1. Ensure that the system is in one of the following states:
 - Ready
 - Warming Up
 - Diagnostic
 - Mechanics Off
 - Check Status
2. Pull out the waste and water drawer to access the water trap
3. Unscrew and remove the water trap. The liquid in the water trap is waste condensation
4. Remove any liquid from the water trap
5. Ensure that the float ball moves freely in the water trap. Ensure that the inner seal of the lid is not missing or damaged
6. Tightly install the water trap
7. Close the waste and water drawer



Monthly Maintenance


Cleaning the Primary Reagent Probes

1. Turn the system mechanics off, and then use the key to unlock the top cover. Lift open the cover and the reagent probes cover
2. To clean the primary reagent probes, push the reagent probe assemblies toward the incubation ring by pushing the motors
3. Wipe each probe with lint-free tissue with cleaning solution, followed by thoroughly rinsing each probe with lint-free tissue with reagent water
4. Reinstall the reagent probes cover, close the top cover, and turn the system mechanics back on



Cleaning the Ancillary Probe and Aspirate Probes

1. Turn the system mechanics off, and then use the key to unlock the top cover. Lift open the top cover
2. To clean the ancillary reagent probe, move the inprocess queue to the right and then pull the probe toward you by holding the top of the probe assembly
3. Clean the probe using a gentle downward motion with lint-free tissue with cleaning solution, followed by thoroughly rinsing the probe with lint-free tissue with reagent water
4. Clean the exterior of the aspirate probes with a lint-free tissue dampened with cleaning solution, followed by rinsing the exterior of the probes with DI water and a lint-free tissue
5. Close the top cover and turn the system mechanics back on



Cleaning the Air Filter


1. Remove the air filter by lifting the cover up and then pulling the cover out of the bracket
2. Remove the old filter and then install the spare, dry air filter
3. Put the cover back onto the system by pushing in and down at the same time
4. Rinse the old air filter with water and then allow it to thoroughly dry



As Needed Maintenance


Cleaning the Reagent Probe Shutter

1. Turn the system mechanics off, and then use the key to unlock the top cover. Lift open the cover and the reagent probes cover
2. Remove the primary reagent packs and refrigerate them
3. Push the reagent probe assemblies toward the incubation ring by pushing the motors
4. Remove the reagent probe shutter by lifting the shutter drive arm up and off the shutter drive and sliding the shutter to the left. Lift the shutter out of the system
5. Wipe the shutter with lint-free tissue with cleaning solution, then rinse with water and thoroughly dry
6. Wipe the surface under the shutter with lint-free tissue with cleaning solution, then rinse with water and thoroughly dry
7. Install the shutter by sliding it to the right until the screw slots are under the screws. Place the shutter arm down on the shutter drive
8. Replace the reagent probes cover, close the top cover, and turn the system mechanics on



Cleaning the Probe Rinse Stations

1. Turn the system mechanics off, and then use the key to unlock the top cover. Lift open the cover and the reagent probes cover
2. Push the reagent, ancillary, and sample probe assemblies toward the incubation ring by pushing the motors
3. Use a transfer pipette to add ~1mL of cleaning solution to each rinse station, then swirl with small maintenance brush. Remove the cleaning solution with a transfer pipette
4. Clean the rinse station with lint-free tissue with cleaning solution, and then rinse the rinse station with lint-free tissue with reagent water
5. Reinstall the reagent probes cover, close the top cover, and turn the system mechanics back on



Cleaning the Primary Reagent Compartment

1. Turn the system mechanics off, and then use the key to unlock the top cover. Lift open the cover and the reagent probes cover
2. Remove the primary reagent packs from the system and refrigerate them
3. Wipe excess water from the primary reagent compartment, and then wipe the surfaces of the compartment with lint-free tissue with water
4. If reagent has spilled onto the reagent pack holders, wipe them with lint-free tissue with water
5. Dry the primary reagent compartment
6. Reinstall the reagent probes cover, close the top cover, and turn the system mechanics back on



Cleaning the Cuvette Waste Chute

1. Open the cuvette waste area door, and remove empty sample tip trays and covers from the tip tray bin
2. Lift the tip tray bin up and then out of the system.
3. Remove the cuvette waste chute by releasing it from the magnetic fitting. Slightly rotate the cuvette waste chute, then pull it out and down
4. Clean the cuvette waste chute with the large maintenance brush and cleaning solution, and then rinse the chute with DI water
5. Dry the inside of the chute with paper towels
6. Install the cuvette waste chute by place the locator pins at the bottom of the chute in the reservoir. Lift the top of the chute until it attaches to the magnetic fitting
7. Install the sample tip tray bin. Install the cuvette waste bin and close the door



Cleaning the Sample Tip Waste Tube

1. Turn the system mechanics off, and then use the key to unlock the top cover. Lift open the cover, and push the sample and ancillary probes back toward the incubation ring
2. Loosen the knurled knob on the sample tip remover, and lift the sample tip remover and sample tip waste tube out of the system
3. Unscrew the sample tip waste tube from the sample tip remover
4. Soak the sample tip remover and sample tip waste tube in cleaning solution for 5 minutes
5. Clean the inside of the sample tip waste tube with cleaning solution and the large maintenance brush
6. Rinse the sample tip remover and sample tip waste tube with water, and then dry thoroughly
7. Connect the sample tip remover to the sample tip waste tube
8. With the slot on the sample tip remover facing away from you, place the bottom of the sample tip waste tube into the opening below the sample probe. Align the 2 holes on the right side with the locator pins on the mounting block. Turn the knurled knob until the sample tip remover is secure
9. Close the top cover, and turn the system mechanics back on



Cleaning the Sample Tip Waste Chute

1. Turn the system mechanics off, and then open the sample tip waste area door and remove the sample tip waste bin
2. Pull out the drawer and lift the lid to the cleaning solution. Ensure the lid remains in the upright position. Wipe the cleaning solution tubing with fresh lint-free tissue to remove excess cleaning solution. Remove the cleaning solution bottle
3. Remove the sample tip waste chute
4. Soak the sample tip waste chute in cleaning solution for 5 minutes
5. While the sample tip waste chute is soaking, wipe the sample tip waste reservoir with lint-free tissue with cleaning solution, and then wipe the reservoir with lint-free tissue with water. Dry the sample tip waste reservoir with lint-free tissue
6. When the sample tip waste chute has finished soaking, clean the chute with the large maintenance brush, and then rinse with water and thoroughly dry the chute
7. Install the sample tip waste chute onto the system with the bottom of the chute in the reservoir
8. Install the cleaning solution bottle and the sample tip waste bin. Close the sample tip waste area door, and turn the system mechanics back on



Cleaning the Covers, Keyboard, and Mouse

1. Carefully clean the following components with lint-free tissue with cleaning solution:
 - Sample entry queue
 - Exterior covers
 - Sample exit queue
 - Ancillary entry
 - Stat entry
 - Keyboard
 - Mouse
 - System fluids tray
 - Tray under the system fluids reservoirs
2. Rinse the components with lint-free tissue with water
3. Dry the components thoroughly

Cleaning the Sample Barcode Scanner


1. Turn the system mechanics off
2. Access the scanner by reaching your hand over and behind the housing
3. Carefully wipe the scanner with an alcohol prep pad
4. Turn the system mechanics back on



Troubleshooting


Clearing Cuvette Obstructions

1. Access the Operator Event Log by selecting **Events** at the Command Bar. Any of these events may lead to checking the orientation and vertical chutes for cuvette jams:
 - Cuvette pusher is offline
 - Cuvette bin is empty (provided the cuvette bin isn't actually empty)
2. Turn the system mechanics off, and then use the key to unlock the top cover. Lift open the cover
3. Observe the cuvette loading bin. Verify that it is not empty
4. Remove any obstructions from the orientation chute
5. Check the vertical chute:
 - Unscrew the vertical chute captive screw and remove the window
 - Remove any obstructions
 - Reinstall the vertical chute window
6. Turn the system mechanics on



Clearing Preheater Jams

1. Access the Operator Event Log by selecting **Events** at the Command Bar. Any of these events may lead to checking the preheater for cuvette jams:
 - Cuvette pusher is offline
 - Ring Loader is offline
2. Turn the system mechanics off, and then use the key to unlock the top cover. Lift open the cover and the reagent probes cover
3. Loosen the thumbscrew on the preheater cover and remove the preheater cover
4. Remove any jammed cuvettes. Ensure that there are 14 cuvettes in the preheater and that the cuvettes do not overlap
5. Replace the preheater cover and tighten the thumbscrew
6. Replace the reagent probes cover, close the top cover, and turn the system mechanics on



Resolving Reagent Probe Errors

1. Access the Operator Event Log by selecting **Events** at the Command Bar. Any of these events may lead to troubleshooting reagent errors:
 - Reagent Probe Volume Check Failure
 - Reagent Probe Vertical Move Error
2. Check the reagent pack for proper placement and volume. The reagent pack that is affected will be displayed in the event message
3. If necessary, continue by removing the pack and checking the probe puncture point. If the point is near the edge of the pack opening, you may need to replace the reagent probe
4. If the reagent pack placement and puncture point are not the issue, perform a dispense test in Diagnostic Tools. If the stream of water is not straight and steady, you may need to replace the reagent probe




Replacing Reagent Probes

1. Turn the system mechanics off, and then use the key to unlock the top cover. Lift open the cover and the reagent probes cover
2. Lift off the reagent probe tubing
3. Use pliers to loosen the reagent probe sleeve
4. Turn the sleeve clockwise, then pull the probe and the sleeve down into the rinse station. Remove the probe and sleeve from the system
5. Place a new reagent probe in the sleeve
6. Lift the probe up into the heater coil and turn it counter clockwise. Tighten the sleeve one quarter turn with pliers
7. Connect the tubing to the top of the probe
8. Replace the reagent probes cover, close the top cover, and turn the system mechanics on



Replacing Aspirate Probes

1. Turn the system mechanics off, and then use the key to unlock the top cover. Lift open the cover
2. Move the inprocess queue to the right
3. Locate the probe that will be replaced and lower the adjacent probe to improve access
4. Disconnect the tubing on the probe being replaced. If the grounded aspirate probe is being replaced, carefully remove the grounding wire from the probe
5. Push the probe pin down and turn counter clockwise to release the spring-loaded probe
6. Remove and discard the aspirate probe
7. Push the new probe down into the probe guide and turn the pin clockwise to secure the spring-loaded probe
8. Connect the probe tubing. If the grounded aspirate probe is being replaced, connect the grounding wire
9. Raise the adjacent probe back up
10. Close the top cover, and turn the system mechanics on



System Restart

Restarting the System

1. At the Status Bar on the bottom of the screen, select System State
2. Select Turn System Off
3. Wait while the system powers down
4. Turn the main power switch off, wait 30 seconds, and then turn the main power switch back on
5. When the software has restarted, sign into the system
6. Turn the mechanics on
7. Wait until the system state is ready, and then return to normal operation



Note: This document is for supplemental use only, and not meant to be used in place of primary technical materials.

© 2014 Siemens Healthcare Diagnostics Inc. All rights reserved.

This guide and the software described within, are copyrighted. No part of this may be copied, reproduced, translated, or reduced to any electronic medium or machine-readable form without the prior written consent of Siemens Healthcare Diagnostics.

All trademarks are the property of their respective owners.

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91052 Erlangen
Germany

Global Division

Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591-5005
USA
www.siemens.com/diagnostics