

Increasing the Sensitivity of MRM Measurements for Tryptic Peptides Using Xevo TQ-S

GOAL

To evaluate the performance of the Xevo™ TQ-S tandem quadrupole mass spectrometer incorporating StepWave™ Technology for tryptic peptide quantification in nanoflow LC-MRM mode.

BACKGROUND

Tandem quadrupole-based Multiple Reaction Monitoring (MRM) mass spectrometry is the most sensitive and robust method for the quantification of target analytes, particularly those in complex biological matrices. However, for the quantification of putative protein biomarkers, sensitivity and specificity are paramount because the target analytes are often present at low concentrations in the relevant tissue or biofluid sample. For this reason, increasing the sensitivity of the MRM assay is crucial for providing a more robust analytical method, which requires less biological material or sample pre-fractionation.

THE SOLUTION

A standard protein digest containing Enolase (P00924), ADH (P00330), BSA (P02769), and Phosphorylase B (P00489) (Part no. 186002865) was dissolved to 1 fmol/μL (H₂O + 0.1% HCO₂H). A Waters® nanoACQUITY UPLC® System was used to separate peptides on a 1.7 μM ACQUITY UPLC® BEH Column (130 C18 75μM x 200 mm) utilizing a gradient from 1% to 40%

Xevo TQ-S provides a highly sensitive and easily automated platform for high-throughput verification of candidate biomarkers.




Figure 1. StepWave ion optics for enhanced sensitivity.

acetonitrile over 90 min at a flow rate of 300 nL/min. The UPLC® eluent was passed directly to the NanoFlow™ Ion Source of the mass spectrometer, which incorporates the StepWave device, as shown in Figure 1.

Two proteotypic peptides from each protein were selected, with the two most intense MRM transitions per peptide chosen for MRM. A total of four MRM transitions per protein were acquired. MRM acquisitions were performed on the Xevo TQ and Xevo TQ-S mass spectrometers. MRM chromatograms for peptide VLYPNENFDGK; m/z 722 to 856, are shown in Figure 2. It can clearly be seen that the signal intensity was significantly increased on the Xevo TQ-S. The absolute sensitivity gain for the MRM transition from 722 to 1067 with the Xevo TQ-S was 73x, as shown in Figure 3, which equates to an eight time improvement in signal-to-noise over the Xevo TQ.

Protein	Peptide Sequence	Peptide transition and ion type	Signal intensity improvement
ADH (P00330)	IGDYAGIK	419 → 666 (y ₆)	11
		419 → 723 (y ₇)	13
	VVGLSTLPEIYEK	724 → 778 (y ₆)	74
		724 → 1080 (y ₆)	62
BSA (P02769)	AEFVEVTK	462 → 476 (y ₂)	52
		462 → 722 (y ₆)	29
	LVNELTEFAK	582 → 595 (y ₆)	37
		582 → 951 (y ₆)	47
Phos B (P00489)	VIFLENYR	527 → 581 (y ₂)	15
		527 → 841 (y ₆)	42
	VLYPNENFDGK	722 → 856 (y ₇)	47
		722 → 1067 (y ₆)	73
Enolase (P00924)	NVNDVIAPAFVK	644 → 745 (y ₇)	47
		644 → 1074 (y ₁₀)	59
	VNQIGTLESISK	645 → 834 (y ₆)	51
		645 → 947 (y ₆)	35

Table 1.

The sensitivity gain for the eight peptides monitored using 16 MRM channels is presented in Table 1.

SUMMARY

Xevo TQ-S shows excellent sensitivity for tryptic peptides at NanoFlow flow rates. Significant improvements in absolute signal intensity and in signal-to-noise ratios for the tryptic peptides were observed in comparison to the Xevo TQ MS. In combination with Verify^E bioinformatics, the nanoACQUITY UPLC System coupled with the Xevo TQ-S provides a highly sensitive and easily automated platform for high-throughput verification of candidate biomarkers.




Figure 2. MRM chromatograms for peptide VLYPNENFDGK; m/z 722 > 856 from Xevo TQ-S (lower panel) and Xevo TQ (top panel).




Figure 3. MRM chromatograms for peptide VLYPNENFDGK, from Glycogen Phosphorylase B; m/z 722 > 1067, from Xevo TQ-S (red) and Xevo TQ (blue) vertical axes are linked to show relative intensities.

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

Waters, UPLC, ACQUITY UPLC, and nanoACQUITY UPLC are registered trademarks of Waters Corporation. The Science of What's Possible, Xevo, StepWave, and NanoFlow are trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2010 Waters Corporation. Produced in the U.S.A.
May 2010 720003491 EN AO-PDF

Waters Corporation
34 Maple Street
Milford, MA 01757 U.S.A.
T: 1 508 478 2000
F: 1 508 872 1990
www.waters.com

