



Flavors and aromas

Analysis of a beefsteak plant

Application Note

Food Testing & Agriculture

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography with an Agilent CP-Wax 52 CB column separates ten components in a sample of a beefsteak plant (*Perilla* sp.) in 40 minutes.



Agilent Technologies

Conditions

Technique : GC-TCT (Off-Line Purging on Tenax TA TCT Tube)
Column : Agilent CP-Wax 52 CB, 0.25 mm x 50 m fused silica
WCOT CP-Wax 52 CB (df = 0.2 µm)
(Part no. CP7723)
Temperature : 80 °C → 220 °C, 5 °C/min
Injector : TCT
Cold trap : CP-Sil 8 CB
0.53 mm; df = 5 µm
Purge time : 10 min
Purge flow : 50 mL/min
Purge temp. : 25 °C
Cryofocussing : -130 °C
Precool time : 3 min
Desorb time : 15 min
Desorb flow : 10 mL/min
Desorption oven temp. : 220 °C
Injection temp. : 210 °C
Injection time : 1 min
GC-injection block temp. : 250 °C


Detector : FID
T = 250 °C

Sample Size : 20 g
Sample Temperature : 40 °C
Concentration Range : ppb/ppm

Courtesy : GL Sciences, Japan

Peak identification

1. acetone
2. α-pinene
3. toluene
4. β-pinene
5. D-limonene
6. linalool
7. β-caryophyllene
8. α-caryophyllene
9. α-farnesene
10. perillaldehyde



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01031



Agilent Technologies