

## Electrophoresis




Ready-to-Run Buffers and Solutions





Bio-Rad is a premier provider of buffers and premixed reagents for life science research. We offer a variety of different products for all your protein and nucleic acid experiments. Whether you need powdered reagents or premixed solutions, Bio-Rad reagents meet the highest quality standards to ensure consistency and reliability in your experiments.



### Electrophoresis Buffers

With premixed electrophoresis running buffers, you standardize your electrophoresis runs and save on preparation time, while avoiding mistakes in buffer concentration. Bio-Rad buffers are made with high-purity water and our own pure reagents, and are 0.4  $\mu\text{m}$  filtered, ensuring the highest quality. Premixed buffers are available for a variety of protein and nucleic acid electrophoresis protocols. Our 5 L boxes offer tremendous economical and convenience advantages. They are compact and stackable to save benchspace, and are designed with an easy-pour spout.



### Blot Processing Buffers

The processing of blots for protein and nucleic acid detection is now even simpler with a variety of premixed wash buffers and blocking solutions.

- Premixed blocking buffers, available as TBS/casein and PBS/casein, take the time and effort out of solubilizing casein
- Premixed wash buffers in TBS, PBS, and SSC reduce the number of stock solutions to prepare
- 10% Tween 20 makes pipetting accurate and simple

### Protein Electrophoresis Buffers

Buffer	1x Formulation	Applications
10x Tris/glycine/SDS	25 mM Tris, 192 mM glycine, 0.1% SDS, pH 8.3	General SDS-PAGE
10x Tris/glycine	25 mM Tris, 192 mM glycine, pH 8.3	Native PAGE
10x Tris/Tricine/SDS	100 mM Tris, 100 mM Tricine, 0.1% SDS, pH 8.3	Peptide SDS-PAGE
10x IEF anode buffer	7 mM phosphoric acid	Analytical isoelectric focusing
10x IEF cathode buffer	20 mM lysine, 20 mM arginine	Analytical isoelectric focusing
10x zymogram renaturation buffer	2.5% Triton X-100	Protease analysis; renatures enzymes after electrophoresis
10x zymogram development buffer	50 mM Tris-HCl, pH 7.5, 200 mM NaCl, 5 mM $\text{CaCl}_2$ , 0.02% Brij 35	Protease analysis; activates enzymes after electrophoresis

### Nucleic Acid Electrophoresis Buffers

10x TBE	89 mM Tris, 89 mM boric acid, 2 mM EDTA, pH 8.3	Nucleic acid electrophoresis/sequencing; polyacrylamide or agarose gels
10x TBE extended range	130 mM Tris, 45 mM boric acid, 2.5 mM EDTA	Nucleic acid electrophoresis/sequencing; polyacrylamide or agarose gels; extends the buffer capacity for longer DNA sequencing runs
50x TAE	40 mM Tris, 20 mM acetic acid, 1 mM EDTA, pH 8.0	Nucleic acid electrophoresis; polyacrylamide or agarose gels

### Blot Processing Buffers

Buffer	1x Formulation	Applications
10x PBS	10 mM sodium phosphate, 150 mM NaCl, pH 7.4	Western blotting wash solution
10x TBS	20 mM Tris, 500 mM NaCl, pH 7.4	Western blotting wash solution, recommended when using alkaline phosphatase
1x PBS/1% casein	10 mM sodium phosphate, 150 mM NaCl, pH 7.4, containing 1% (w/v) casein	Western blotting blocking buffer; casein blockers recommended for all applications, including those with biotin-avidin complexes
1x TBS/1% casein	20 mM Tris, 500 mM NaCl, pH 7.4, containing 1% (w/v) casein	Western blotting blocking buffer; casein blockers recommended for all applications, including those with biotin-avidin complexes
20x SSC	150 mM sodium chloride, 15 mM sodium citrate, pH 7.0	Northern and Southern blotting prehybridization and hybridization solutions
Tween 20	10% w/v Tween 20 or 100% Tween 20	Blocking and wash buffer component



### Sample Loading Buffers

Premixed loading buffers remove variables that cause lane-to-lane running anomalies, and since no preparation is required, you save valuable time as well. Bio-Rad premixed sample buffers are available for numerous applications, including native PAGE, SDS-PAGE, peptide analysis, analytical IEF, nucleic acid sample preparation (denaturing and nondenaturing), and zymogram gel sample preparation.



### Blot Transfer Buffers

The transfer buffer must facilitate both effective elution from the gel matrix and effective binding of the protein or nucleic acid to the membrane. Determine your choice of buffer by the type of gel or membrane and the physical characteristics of the molecules of interest.

#### Sample Loading Buffers

Buffer	Formulation	Applications
Laemmli sample buffer	62.5 mM Tris-HCl, pH 6.8, 2% SDS, 25% glycerol, 0.01% Bromophenol Blue	SDS-PAGE
Native sample buffer	62.5 mM Tris-HCl, pH 6.8, 40% glycerol, 0.01% Bromophenol Blue	PAGE
Tricine sample buffer	200 mM Tris-HCl, pH 6.8, 2% SDS, 40% glycerol, 0.04% Coomassie G-250	Peptide analysis, small protein SDS-PAGE
IEF sample buffer	50% glycerol	Isoelectric focusing
Zymogram sample buffer	62.5 mM Tris-HCl, pH 6.8, 25% glycerol, 4% SDS, 0.01% Bromophenol Blue	Protease analysis
Nucleic acid sample buffer (5x)	50 mM Tris-HCl pH 8.0, 25% glycerol, 5 mM EDTA, 0.2% Bromophenol Blue, 0.2% Xylene Cyanol FF	Nondenaturing dsDNA
TBE-urea sample buffer	89 mM Tris-HCl, pH 8.0, 89 mM boric acid, 2 mM EDTA, 7 M urea, 12% Ficoll, 0.01% Bromophenol Blue, 0.02% Xylene Cyanol FF	Denaturing ssDNA, RNA

#### Western Blotting

Buffer	1x Formulation	Applications
10x Tris/glycine	25 mM Tris, 192 mM glycine, pH 8.3	<b>SDS-PAGE gels (tank or semi-dry blotting):</b> Add 20% methanol to remove SDS from the protein and improve its affinity for nitrocellulose <b>Native PAGE gels (tank blotting):</b> For acidic and neutral proteins, use Tris/glycine buffer without methanol
10x Tris/CAPS	60 mM Tris, 40 mM CAPS	<b>SDS-PAGE (semi-dry blotting only):</b> Discontinuous buffer system increases transfer efficiency; to Tris/CAPS buffer add 15% methanol for the anode buffer and 0.1% SDS for the cathode buffer

#### Southern/Northern Blotting

50x TAE	40 mM Tris, 20 mM acetic acid, 1 mM EDTA, pH 8.0	Tank blotting of polyacrylamide gels
10x TBE	89 mM Tris, 89 mM boric acid, 2 mM EDTA, pH 8.3	Tank blotting or semi-dry blotting of polyacrylamide gels
20x SSC	150 mM sodium chloride, 15 mM sodium citrate, pH 7.0	Capillary transfer of agarose gels

If you don't find exactly what you need, simply contact your local Bio-Rad representative and inquire about custom-made buffers.



### SDS Solutions

Detergents are employed in electrophoresis when it is necessary to disrupt protein-lipid or protein-protein interactions. SDS is the most common detergent used in PAGE analysis because most proteins are readily soluble in it. Bio-Rad SDS solutions are highly purified — an important feature, since impurities in SDS have unpredictable effects on electrophoretic mobilities.

### Gel Casting Solutions

#### Tris Buffers and Acrylamide Solutions for Gel Casting

Bio-Rad offers a variety of prepared solutions for casting polyacrylamide gels. Tris solutions are formulated into working concentrations for preparing the stacking and resolving portions of native or SDS-PAGE gels, according to Laemmli or Ornstein-Davis discontinuous buffer systems. Acrylamide solutions are provided ready to use and come with instructions. High-purity reagents and carefully controlled manufacturing conditions allow acrylamide solutions to be stable for 1 year at 4°C.



### Electrophoresis Buffer Reagents

In case you would like to prepare it all yourself, we offer a complete line of reagents. Our classic electrophoresis powder reagents are the ultimate in high quality.

### SDS Solutions

Solution	Formulation	Applications
10% SDS solution	10% (w/v) sodium dodecyl sulfate	SDS-PAGE: for preparing sample, gel, and running buffers
20% SDS solution	20% (w/v) sodium dodecyl sulfate	Northern and Southern hybridization buffer component

### Gel Casting Buffers

Solution	Applications
1.5 M Tris-HCl, pH 8.8	Resolving gel preparation
0.5 M Tris-HCl, pH 6.8	Stacking gel preparation
Acrylamide Solutions	
19:1 Acrylamide/Bis	DNA sequencing
29:1 Acrylamide/Bis	Protein separation
37.5:1 Acrylamide/Bis	Protein separation



## Ordering Information

Catalog # Description

### Electrophoresis Running Buffers

161-0732 10x Tris/Glycine/SDS, 1 L  
161-0772 10x Tris/Glycine/SDS, 5 L cube  
161-0734 10x Tris/Glycine, 1 L  
161-0771 10x Tris/Glycine, 5 L cube  
161-0744 10x Tris/Tricine/SDS, 1 L  
161-0761 10x IEF Anode Buffer, 250 ml  
161-0762 10x IEF Cathode Buffer, 250 ml  
161-0765 10x Zymogram Renaturation Buffer, 125 ml  
161-0766 10x Zymogram Development Buffer, 125 ml  
161-0733 10x Tris/Boric Acid/EDTA (TBE), 1 L  
161-0770 10x Tris/Boric Acid/EDTA (TBE), 5 L cube  
161-0741 10x TBE Extended Range, 1 L  
161-0743 50x Tris/Acetic Acid/EDTA (TAE), 1 L  
161-0773 50x Tris/Acetic Acid/EDTA (TAE), 5 L cube

### Blot Processing Buffers

170-6435 10x TBS, 1 L  
161-0780 10x PBS, 1 L  
161-0783 1x PBS/1% Casein, 1 L  
161-0782 1x TBS/1% Casein, 1 L  
161-0774 20x SSC, 1 L  
161-0775 20x SSC, 5 L cube  
161-0781 10% Tween 20, 1 L  
170-6531 Tween 20, 100 ml

### Sample Loading Buffers

161-0737 Laemmli Sample Buffer, 30 ml  
161-0738 Native Sample Buffer, 30 ml  
161-0739 Tricine Sample Buffer, 30 ml  
161-0763 IEF Sample Buffer, 30 ml  
161-0764 Zymogram Sample Buffer, 30 ml  
161-0767 Nucleic Acid Sample Buffer, 5x, 10 ml  
161-0768 TBE-Urea Sample Buffer, 30 ml

### Blotting Transfer Buffers

161-0734 10x Tris/Glycine, 1 L  
161-0771 10x Tris/Glycine, 5 L cube  
161-0778 10x Tris/CAPS, 1 L  
161-0743 50x TAE, 1 L  
161-0773 50x TAE, 5 L cube  
161-0733 10x TBE, 1 L  
161-0770 10x TBE, 5 L cube  
161-0774 20x SSC, 1 L  
161-0775 20x SSC, 5 L cube

Catalog # Description

### SDS Solutions

161-0416 SDS Solution, 10%, 250 ml  
161-0418 SDS Solution, 20%, 1,000 ml

### Buffer Reagents

161-0716 Tris, 500 g  
161-0719 Tris, 1 kg  
161-0729 EDTA, 500 g  
161-0717 Glycine, 250 g  
161-0718 Glycine, 1 kg  
161-0724 Glycine, 2 kg  
161-0713 Tricine, 500 g  
161-0730 Urea, 250 g  
161-0731 Urea, 1 kg  
161-0610 Dithiothreitol, 1 g  
161-0611 Dithiothreitol, 5 g  
161-0710 2-Mercaptoethanol, 25 ml  
161-0301 SDS, 100 g  
161-0302 SDS, 1 kg

### Gel Casting Solutions

161-0798 1.5 M Tris-HCl, pH 8.8, 1 L  
161-0799 0.5 M Tris-HCl, pH 6.8, 1 L  
161-0154 30% Acrylamide/Bis Solution 19:1, 500 ml  
161-0156 30% Acrylamide/Bis Solution 29:1, 500 ml  
161-0158 30% Acrylamide/Bis Solution 37.5:1, 500 ml  
161-0144 40% Acrylamide/Bis Solution 19:1, 500 ml  
161-0146 40% Acrylamide/Bis Solution 29:1, 500 ml  
161-0148 40% Acrylamide/Bis Solution 37.5:1, 500 ml

Brij and Tween are trademarks of ICI Americas Inc. Ficoll is a trademark of GE Healthcare group companies. Triton is a trademark of Union Carbide. Coomassie is a trademark of BASF Aktiengesellschaft.

**BIO-RAD**

**Bio-Rad  
Laboratories, Inc.**

Life Science  
Group

Web site [www.bio-rad.com](http://www.bio-rad.com) USA 800 4BIORAD Australia 61 02 9914 2800 Austria 01 877 89 01 Belgium 09 385 55 11 Brazil 55 21 3237 9400  
Canada 905 364 3435 China 86 21 6426 0808 Czech Republic 420 241 430 532 Denmark 44 52 10 00 Finland 09 804 22 00 France 01 47 95 69 65  
Germany 089 318 84 0 Greece 30 210 777 4396 Hong Kong 852 2789 3300 Hungary 36 1 455 8800 India 91 124 4029300 Israel 03 963 6050  
Italy 39 02 216091 Japan 03 6361 7000 Korea 82 2 3473 4460 Mexico 52 555 488 7670 The Netherlands 0318 540666 New Zealand 0508 805 500  
Norway 23 38 41 30 Poland 48 22 331 99 99 Portugal 351 21 472 7700 Russia 7 495 721 14 04 Singapore 65 6415 3188 South Africa 27 861 246 723  
Spain 34 91 590 5200 Sweden 08 555 12700 Switzerland 061 717 95 55 Taiwan 886 2 2578 7189 United Kingdom 020 8328 2000