



CHAN ZUCKERBERG
BIOHUB

UCSF

University of California
San Francisco

Metagenomic detection of antimicrobial resistant respiratory pathogens

Chaz Langelier, MD, PhD

Assistant Professor, Division of Infectious Diseases


Assoc. Medical Director, Hospital Epidemiology & Infection Prevention

University of California San Francisco


Chan Zuckerberg Biohub

Lower respiratory tract infections (LRTI) are a leading cause of global mortality


Infectious Disease Mortality in the USA




Leading Causes of Death Worldwide






A need for better LRTI diagnostics



Mortality from Antimicrobial Resistant Infections



Current paradigm of empiric treatment for LRTI contributes to the AMR problem

- **Antibiotic overtreatment (40-67% cases)** 
 - Selection for AMR pathogens
 - Microbiome disruption: *C. difficile* infection 
 - Adverse drug reactions/allergies 
 - Increased morbidity and mortality
- **Antibiotic undertreatment**
 - Infection progression and/or transmission

Metagenomics: a key tool to combat AMR


- **Early detection**

- Effective treatment of resistant microbes to avoid treatment failure
- Targeted therapy to improve antimicrobial stewardship

- **Surveillance**

- Global transmission patterns
- Regional distribution of resistant pathogens
- Emergence of novel resistance mechanisms

mNGS for **AMR** pathogen detection



mNGS vs. gold standard of culture-based phenotypic antimicrobial susceptibility testing

1. mNGS alone (RNA-Seq + DNA-Seq)
2. mNGS + Cas9 targeted enrichment
3. Nanopore + Cas9 targeted enrichment
4. Antimicrobial stewardship modeling

mNGS for **AMR** pathogen detection


Phenotypic
Testing:

 Resistant

 Susceptible


		Pen/Amp	Vancomycin	Methicillin	Clindamycin	TMP/SMX	Ceftriaxone
225	<i>Staphylococcus aureus</i>	Resistant	Susceptible	Resistant	Resistant	Susceptible	n/a
288	<i>Staphylococcus aureus</i>	Resistant	Susceptible	Resistant	Resistant	Susceptible	n/a
409	<i>Staphylococcus aureus</i>	Resistant	Susceptible	Resistant	Resistant	Susceptible	n/a
212	<i>Staphylococcus aureus</i>	Resistant	Susceptible	Susceptible	Resistant	Susceptible	n/a
230	<i>Staphylococcus aureus</i>	Resistant	Susceptible	Resistant	Resistant	Susceptible	n/a
239	<i>Staphylococcus aureus</i>	Resistant	Susceptible	Resistant	Resistant	Susceptible	n/a
289	<i>Staphylococcus aureus</i>	Resistant	Susceptible	Susceptible	Susceptible	Susceptible	n/a
304	<i>Staphylococcus aureus</i>	Resistant	Susceptible	Susceptible	Susceptible	Susceptible	n/a
311	<i>Staphylococcus aureus</i>	Resistant	Susceptible	Susceptible	Susceptible	Susceptible	n/a
314	<i>Staphylococcus aureus</i>	Resistant	Susceptible	Susceptible	Susceptible	Susceptible	n/a
386	<i>Staphylococcus aureus</i>	Resistant	Susceptible	Susceptible	Susceptible	Susceptible	n/a
407	<i>Staphylococcus aureus</i>	Resistant	Susceptible	Susceptible	Susceptible	Susceptible	n/a
277	<i>Enterococcus faecium</i>	Resistant	Resistant	n/a	Resistant	n/a	n/a
311	<i>Streptococcus pneumoniae</i>	Susceptible	Susceptible	n/a	n/a	n/a	Susceptible
334	<i>Streptococcus pneumoniae</i>	Susceptible	Susceptible	n/a	n/a	n/a	Susceptible
337	<i>Streptococcus pneumoniae</i>	Resistant	Susceptible	n/a	n/a	n/a	Resistant
399	<i>Streptococcus pneumoniae</i>	Susceptible	Susceptible	n/a	n/a	n/a	Susceptible
463	<i>Streptococcus pneumoniae</i>	Susceptible	Susceptible	n/a	n/a	n/a	Susceptible
213	<i>Escherichia coli</i>	n/a	n/a	n/a	n/a	Susceptible	Susceptible
213	<i>Klebsiella pneumoniae</i>	n/a	n/a	n/a	n/a	Susceptible	Susceptible
232	<i>Klebsiella pneumoniae</i>	Resistant	n/a	n/a	n/a	Susceptible	Resistant
251	<i>Klebsiella pneumoniae</i>	n/a	n/a	n/a	n/a	Susceptible	Susceptible
239	<i>Pseudomonas aeruginosa</i>	Resistant	n/a	n/a	n/a	Susceptible	Resistant
268	<i>Pseudomonas aeruginosa</i>	Resistant	n/a	n/a	n/a	Susceptible	Resistant
268	<i>Stenotrophomonas</i>	Resistant	n/a	n/a	n/a	Susceptible	Resistant
257	<i>Enterobacter aerogenes</i>	Resistant	n/a	n/a	n/a	Susceptible	Resistant
298	<i>Enterobacter cloacae</i>	Resistant	n/a	n/a	n/a	Susceptible	Resistant
343	<i>Burkholderia cepaciae</i>	n/a	n/a	n/a	n/a	Susceptible	Susceptible
350	<i>Serratia marcescans</i>	Resistant	n/a	n/a	n/a	Susceptible	Resistant
382	<i>Serratia marcescans</i>	Resistant	n/a	n/a	n/a	Susceptible	Resistant
485	<i>Haemophilus influenzae</i>	Resistant	n/a	n/a	n/a	Susceptible	Susceptible

mNGS methods for **AMR** pathogen detection



Rules-based AMR model:

- SRST2 + custom AMR gene database



- 5% minimum allele coverage
- SNP detection: CARD RGI tool
- Phenotypic resistance prediction:
 - CARD ontology
 - 2 MD adjudication

mNGS vs culture for AMR prediction

S. aureus only (n = 12)

	Sensitivity	Specificity	Accuracy
DNA-Seq	52%	92%	80%
RNA-Seq	40%	95%	77%
combined	83%	95%	91%

All (n = 31)

	Sensitivity	Specificity	Accuracy
DNA-Seq	52%	92%	80%
RNA-Seq	69%	96%	88%
combined	88%	92%	91%

A 'false-positive' case

- 2/25 56 y/o man admitted for venous thrombosis
- 2/27 New fever, hypotension, hypoxia
 - Rx: aztreonam, vancomycin, azithromycin
- 2/28 Transferred to ICU, respiratory Cx obtained
 - Enrolled in study, sample for mNGS collected
- 3/02 Cx grows *Klebsiella*
 - Vancomycin, azithromycin stopped
- 3/04 Improved, discharged from ICU
- 3/07 new hypoxemic respiratory failure
 - Re-admitted to ICU, Rx vancomycin, TA Cx obtained
- 3/13 TA Cx finalize: *MRSA*

mNGS:
-*K. pneumoniae*

AMR genes:
-*SHV-190*
-*mecA*
-*ermA/C*

Staph. aureus

Susceptibility testing:

Clindamycin	Resistant
Nafcillin	Resistant
Penicillin G	Resistant
Bactrim	Susceptible
Vancomycin	Susceptible

mNGS vs culture for AMR prediction


S. aureus only (n = 12)

	Sensitivity	Specificity	Accuracy
DNA-Seq	52%	92%	80%
RNA-Seq	40%	95%	77%
combined	83%	95%	91%


All (n = 31)

	Sensitivity	Specificity	Accuracy
DNA-Seq	52%	92%	80%
RNA-Seq	69%	96%	88%
combined	88%	92%	91%


Resistance gene sequences are rare, ~ 1 in 10^8




FLASH: Finding Low Abundance Sequences by Hybridization




FLASH CRISPR/Cas9 guide RNA design




FLASH of clinical *Staphylococcus* isolates



FLASH: how many reads are needed?




FLASH limit of detection




Cultured MRSA Isolates

FLASH: Is there an upper limit to guide RNA multiplexing?




FLASH Cas9 Targeted mNGS - Example



FLASH-mNGS on direct clinical samples

Staph aureus LRTI



Tracheal aspirate

mNGS + FLASH vs culture for AMR prediction

S. aureus only (n = 12)

	Sensitivity	Specificity	Accuracy
DNA-Seq	52%	92%	80%
RNA-Seq	40%	95%	77%
combined	83%	95%	91%
+ FLASH	95%	93%	93%


All (n = 31)

	Sensitivity	Specificity	Accuracy
DNA-Seq	52%	92%	80%
RNA-Seq	69%	96%	88%
combined	88%	92%	91%
+ FLASH	97%	89%	92%

Rapid mNGS-FLASH via Nanopore?




Rapid mNGS-FLASH via Nanopore





Wayne Deng PhD; Charles Chiu, MD PhD

Precision Antibiotic Stewardship




- Median time to AST: 93 hours
- Empiric antibiotics: 91% of patients
- Theoretical impact of mNGS:
 - Adjudicated optimal Δ antibiotics based on 48h mNGS turnaround
 - Compared against actual antibiotics received during hospitalization
 - **39% reduction** in unnecessary days of antibiotic therapy


ICU Respiratory Resistome



Epidemiological classification of pneumonia




Epidemiological classification of pneumonia




**Recent guidelines
group together
HCAP and CAP**

AMR gene burden differs based on epidemiological classification of pneumonia



Airway β -lactamase genes in acute trauma patients at time of ICU admission vs 48 hours later



n = 15

Jen Mulliken, MD

Resistome of a hospital room

ESBL and
carbapenemase
genes

β -lactamase genes

OXA-23_Bla	
AIM-1_Bla	←
OXA-211_Bla	
OXA-2_Bla	
OXA-237_Bla	
TEM-1D_Bla	
IND-1_Bla	←
OXA-50_Bla	←
CTX-M-1_Bla	←
BlaZ_Bla	
CfxA_Bla	
CME-1_Bla	←
OCH-2_Bla	←
GOB-16_Bla	←
CEPH-AQU_Bla	
CMY-MOX_Bla	←
OXA-60_Bla	←

Future mNGS AMR Directions

- **Hospital resistome surveillance**
- **Dynamics of microbiome and AMR genes based on daily antibiotic exposure**
 - Children with acute respiratory illnesses
 - Lung transplant recipients
- **mNGS-FLASH for sepsis pathogen diagnosis**

Acknowledgements



University of California
San Francisco

Katrina Kalantar

Jack Kamm

Emily Crawford

Joe DeRisi

Carolyn Calfee

Charles Chiu

Pete Mourani

Saharai Caldera

Paula Serpa

Josh Batson

Jen Mulliken

Aaron McGeever