





1-800-224-6113 support@intermedtn.com www.intermedtn.com

CLEAN & CLEAR. *EVERY TIME.*™

# CLEAN & CLEAR. EVERY TIME.

The Dornier AXIS™ Single-Use Digital Ureteroscope provides every patient a clean scope with a clear view - every time. Eliminate the cost and risk associated with reprocessing. You and your patient deserve an innovative, quality ureteroscope that helps bring quality care to the next level.

The Dornier  $\mathsf{AXIS}^\mathsf{TM}$  - for a clean & clear scope, every time.



**Traditional look** 

ergonomic design

in a lightweight

and feel

275°
Dual Deflection

Excellent column strength at 8.5F OD

Integrated irrigation port



technology for high resolution

Dual LED illumination

2.5X
more imaging pixels\*
High pixel count

**ECLIPSE**<sup>TM</sup>

Zero Tip Basket

for clear and sharp image



| Product Name                  | Part Number | Product Description                                      |
|-------------------------------|-------------|--|
| Disc™                         | BIPDDSC     | Plasma Disc Electrode, 6/box                             |
| 12 Degree Medium Cutting Loop | BIP12CLM    | 12 Degree Medium Cutting Loop Electrode, 6/box           |
| 30 Degree Medium Cutting Loop | BIP30CLM    | 30 Degree Medium Cutting Loop Electrode, 6/box           |
| 12 Degree Large Cutting Loop  | BIP12CLL    | 12 Degree Large Cutting Loop Electrode, 6/box            |
| 30 Degree Large Cutting Loop  | BIP30CLL    | 30 Degree Large Cutting Loop Electrode, 6/box            |
| Needle Electrode              | BIPNEEL     | Needle Electrode, 6/box                                  |
| Medium Bladder Loop           | BIPMBLA     | Medium Bladder Loop Electrode, 6/box                     |
| Bipolar Cable, ESG-400        | BIPCBESG    | Single-use Bipolar Cable, ESG-400 Generator, 3/box       |
| Bipolar Cable, PK SuperPulse  | BIPCBPKS    | Single-use Bipolar Cable, PK SuperPulse Generator, 3/box |

ESG-400 Compatible Cable



PK SuperPulse Compatible Cable





1-800-224-6113
support@intermedtn.com
www.intermedtn.com



## Where Design Meets Functionality



**True 360° table access** with free-standing treatment table. Unrivaled flexibility through C-arm parking position.



High-resolution imaging from unprecedented angles with the flexible c-arm and table movements. *Independent or with sophisticated synchronization.* 





# **Urology Applications**

The unique design of the Nautilus®, makes this urotable the ideal workstation for a wide range of urological procedures such as:

Laparoscopic Procedures Transurethral Procedures

Diagnostic Urology Procedures Contrast Studies Percutaneous Endourology Procedures Micturition
Cystourethrogram
(MCU)

## **Extending New Dimensions**

**Raising new standards** in radiation management through the flexible positioning of x-ray tube under table.



- 1 Free Table Access
  C-arm parking position with 90° swiveling.
- 2 Motorized Free-standing Table Isocentric and non-isocentric Trendelenburg table tilting.
  Fully radiolucent table top made of carbon fiber.
- 3 43 x 43 FPD with SID-hub
  Detector-to-source flexibility
  captures KUB (Kidney-UreterBladder) in one shot.
  - 4 21" Monitors

    Easy access view to monitors
    from any position; interfacing
    with additional equipment
    (e. g. endoscopes, ultrasound).
- 5 80 kW Advanced High Frequency Generator High efficiency and distinguished performance.
- 6 Accessories Shelf
  Three compartments for surgical equipment.



#### **InterMed Resources TN**

9005 Overlook Blvd. Brentwood, TN 37027 +1-800-224-6113







# SINGLE-USE DIGITAL FLEXIBLE URETEROSCOPE

PU3033A




# Precise & Extraordinary

Reduced Outer Diameter Irrigation Up To

19% 147%

Both Compared with PU3022A in 12/14 Fr Sheath In 10/12 Fr sheath, irrigation is 280% compared with PU3022A In 11/13 Fr sheath, irrigation is 179% compared with PU3022A





≥9 Fr ID Sheath **Adaptable** 



Lower **Complication Rate** 



Higher **Stone-Free Rate** 

| Optical System                                 |                     |  |  |
|--|---------------------|--|--|
| Field of View                                  | 90°                 |  |  |
| Direction of View                              | 0°                  |  |  |
| Depth of Field                                 | 3-50 mm             |  |  |
| Insertio                                       | n Portion           |  |  |
| Controllable Portion                           | 270° Up   270° Down |  |  |
| Maximum Insertion Portion W                    | <b>lidth</b> 2.7 mm |  |  |
| Diameter of Insertion Portion                  | 7.5 Fr              |  |  |
| Distal End Diameter                            | 7.5 Fr              |  |  |
| Working Length                                 | 650 mm              |  |  |
| Minimum Instrument Channe                      | Width 3.6 Fr        |  |  |
| Additional Features                            |                     |  |  |
| <b>Custom Buttons</b>                          | PU3033A ●           |  |  |
| Stepless Self-locking Lever (Deflection Lever) | PU3033A ●           |  |  |



PIONEERING TECHNOLOGIES
AND REVOLUTIONARY
THERAPIES FOR THE
ULTIMATE TREATMENT
EXPERIENCE IN UROLOGY.

CONTACT INTERMED RESOURCES TN TODAY TO CONNECT TO YOUR REPRESENTATIVE FOR FURTHER INFORMATION.

Leading Technology. Improving Life.



YOUR ONE UROLOGY SOURCE



1-800-224-6113 support@intermedtn.com www.intermedtn.com

