

Radial Artery Catheterization Kit



ARROW
INTERNATIONAL

Arrow Select Kits

WBH-04020-1

Contents:

- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: 1% Lidocaine HCl Solution
- 1: Swabstick Foil Pack Skin Protectant (non-sterile solution)
- 1: **1.5 mL Applicator 2% CHG and 70% IPA ChloroPrep^{®1} One-Step Frepp[®] Solution (non-sterile solution)**
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1: Suture: 3-0 Braided Silk with Straight Needle
- 1: Local Anesthetic Components
- 4: Gauze Pads: 4" x 4"
- 1: Safety Scalpel: #11
- 1: SharpsAway IIT[™] Locking Disposal Cup
- 1: Wing Clip: Suture

¹ A registered trademark of CareFusion or one of its subsidiaries.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Single use – do not resterilize. Sterilized by ethylene oxide. Keep away from sunlight. Keep dry. Do not use if package has been previously opened or damaged. Not made with natural rubber latex. Fluid path components are non-pyrogenic. Contains medication.

Store between 20 - 25°C (68 - 77°F).



Arrow International LLC
3015 Carrington Mill Blvd., Morrisville, NC 27560 USA
www.teleflex.com

IH-04020-100G Rev. 3 (2021-09)

Radial Artery Catheterization Kit

WBH-04020-1

Arrow Select Kits



ARROW
INTERNATIONAL

20^{Ga.}

4.45^{cm}
catheter
length

.018^{inch}
diameter
spring-wire
guide