





# Knob weights




## STAINLESS STEEL HF12





highly polished

### SINGLE WEIGHT

Nominal value	Tol.	without box	 Plastic boxes with screwed top	 Clean room plastic cases*	 Mahogany wooden box	 Clean room case with insert	<b>DAkks</b> Calibration certificate
+/- mg		Art.no.	Art.no.	Art.no.	Art.no.	Art.no.	8.0000-Art.no.
1 g	0,030	8.MEHO-110	8.MEHD-110	8.MEHP-110	8.MEHM-110		...111
2 g	0,040	8.MEHO-120	8.MEHD-120	8.MEHP-120	8.MEHM-120		...121
5 g	0,050	8.MEHO-130	8.MEHD-130	8.MEHP-130	8.MEHM-130		...131
10 g	0,060	8.MEHO-140	8.MEHD-140	8.MEHP-140	8.MEHM-140		...141
20 g	0,080	8.MEHO-150	8.MEHD-150	8.MEHP-150	8.MEHM-150		...151
50 g	0,10	8.MEHO-160	8.MEHD-160	8.MEHP-160	8.MEHM-160		...161
100 g	0,16	8.MEHO-170	8.MEHD-170	8.MEHP-170	8.MEHM-170	8.MEHR-170	...171
200 g	0,3	8.MEHO-180	8.MEHD-180	8.MEHP-180	8.MEHM-180	8.MEHR-180	...181
500 g	0,8	8.MEHO-190	8.MEHD-190	8.MEHP-190	8.MEHM-190	8.MEHR-190	...191
1 kg	1,6	8.MEHO-210	8.MEHD-210	8.MEHP-210	8.MEHM-210	8.MEHR-210	...211
2 kg	3,0	8.MEHO-220	8.MEHD-220	8.MEHP-220	8.MEHM-220	8.MEHR-220	...221
5 kg	8,0	8.MEHO-230	8.MEHD-230	8.MEHP-230	8.MEHM-230	8.MEHR-230	...231
10 kg	16	8.MEHO-240	8.MEHD-240	8.MEHP-240	8.MEHM-240	8.MEHR-240	...241
20 kg	30	8.MEHO-250			8.MEHM-250	8.MEHR-250	...251
50 kg	80	8.MEHO-260			8.MEHM-260	8.MEHR-260	...261

Nominal value	Tol.				 Service case	<b>DAkks</b> Calibration certificate
+/- mg					Art.-Nr.	8.0000-Art.-Nr.
5 kg	8,0				8.MEHS-230	...231
10 kg	16				8.MEHS-240	...241
20 kg	30				8.MEHS-250	...251
50 kg	80				8.MEHS-260	...261

### SETS

Nominal value	Qty.	Content	 Clean room plastic box	 Mahogany wooden box	 Service case	 Clean room case with insert	<b>DAkks</b> Calibration certificate
			Art.no.	Art.no.	Art.no.	Art.no.	8.0000-Art.no.
1 g – 50 g	8	110 g	8.XTHP-560	8.XTHM-560		8.XTHR-560	...561
1 g – 100 g	9	210 g	8.XTHP-570	8.XTHM-570	8.XTHS-570	8.XTHR-570	...571
1 g – 200 g	11	610 g	8.XTHP-580	8.XTHM-580	8.XTHS-580	8.XTHR-580	...581
1 g – 500 g	12	1.110 g	8.XTHP-590	8.XTHM-590	8.XTHS-590	8.XTHR-590	...591
1 g – 1 kg	13	2.110 g	8.XTHP-610	8.XTHM-610	8.XTHS-610	8.XTHR-610	...611
1 g – 2 kg	15	6.110 g	8.XTHP-620	8.XTHM-620	8.XTHS-620	8.XTHR-620	...621
1 g – 5 kg	16	11.110 g	8.XTHP-630	8.XTHM-630	8.XTHS-630	8.XTHR-630	...631
1 g – 10 kg	17	21.110 g	8.XTHP-640	8.XTHM-640	8.XTHS-640	8.XTHR-640	...641
1 kg – 5 kg	4	10.000 g	8.XTHP-680	8.XTHM-680	8.XTHS-680	8.XTHR-680	...681

\* up from 100 g in a plastic case with foam insert





# Knob weights







**STAINLESS STEEL HF12**

highly polished

SETS with mg- weights (flat polygonal sheets)

Nominal value	Qty.	Disc weights	 Clean room plastic box	 Mahogany wooden box	 Service case	 Clean room case with insert	<b>DAkks</b> Calibration certificate
		Art.no.	Art.no.	Art.no.	Art.no.	8.0000-Art.no.	
1 mg – 5 g	16	6,11 g	8.XFHP-730	8.XFHM-730		8.XFHR-730	...731
1 mg – 50 g	20	111,11 g	8.XFHP-760	8.XFHM-760		8.XFHR-760	...761
1 mg – 100 g	21	211,11 g	8.XFHP-770	8.XFHM-770	8.XFHS-770	8.XFHR-770	...771
1 mg – 200 g	23	611,11 g	8.XFHP-780	8.XFHM-780	8.XFHS-780	8.XFHR-780	...781
1 mg – 500 g	24	1.111,11 g	8.XFHP-790	8.XFHM-790	8.XFHS-790	8.XFHR-790	...791
1 mg – 1 kg	25	2.111,11 g	8.XFHP-810	8.XFHM-810	8.XFHS-810	8.XFHR-810	...811
1 mg – 2 kg	27	6.111,11 g	8.XFHP-820	8.XFHM-820	8.XFHS-820	8.XFHR-820	...821
1 mg – 5 kg	28	11.111,11 g	8.XFHP-830	8.XFHM-830	8.XFHS-830	8.XFHR-830	...831
1 mg – 10 kg	29	21.111,11 g	8.XFHP-840	8.XFHM-840	8.XFHS-840	8.XFHR-840	...841
1 mg – 1 kg (2x)	38	3.778,22 g	8.XFHP-850	8.XFHM-850	8.XFHS-850	8.XFHR-850	...851

SETS Sets with mg- weights (wire shape)

Nominal value	Qty.	Content	 Clean room plastic box	 Mahogany wooden box	 Service case	 Clean room case with insert	<b>DAkks</b> Calibration certificate
		Art.-Nr.	Art.-Nr.	Art.-Nr.	Art.-Nr.	8.0000-Art.-Nr.	
1 mg – 5 g	16	6,11 g	8.XEHP-730	8.XEHM-730		8.XEHR-730	...731
1 mg – 50 g	20	111,11 g	8.XEHP-760	8.XEHM-760		8.XEHR-760	...761
1 mg – 100 g	21	211,11 g	8.XEHP-770	8.XEHM-770	8.XEHS-770	8.XEHR-770	...771
1 mg – 200 g	23	611,11 g	8.XEHP-780	8.XEHM-780	8.XEHS-780	8.XEHR-780	...781
1 mg – 500 g	24	1.111,11 g	8.XEHP-790	8.XEHM-790	8.XEHS-790	8.XEHR-790	...791
1 mg – 1 kg	25	2.111,11 g	8.XEHP-810	8.XEHM-810	8.XEHS-810	8.XEHR-810	...811
1 mg – 2 kg	27	6.111,11 g	8.XEHP-820	8.XEHM-820	8.XEHS-820	8.XEHR-820	...821
1 mg – 5 kg	28	11.111,11 g	8.XEHP-830	8.XEHM-830	8.XEHS-830	8.XEHR-830	...831
1 mg – 10 kg	29	21.111,11 g	8.XEHP-840	8.XEHM-840	8.XEHS-840	8.XEHR-840	...841
1 mg – 1 kg (2x)	38	3.778,22 g	8.XEHP-850	8.XEHM-850	8.XEHS-850	8.XEHR-850	...851

## 3-Part service case

Available for all classes of weights

1mg – 5kg



Price on request

## Trolley

Available for all classes of weights

1mg – 10kg  
1mg – 20kg



Price on request