

**ONDAS DE CHOQUE (4 BARS, 15 Hz.) – EL.80-CH****ONDAS DE CHOQUE (5 BARS, 20 Hz.) – EL.81-CH**

20 Protocolos memorizados

Memoria interna con 200 programas para memorizar

Voltaje: 220-240V/50-60Hz. Bajo pedido 110-127V/50-60Hz.

EL.80-CH: Presión hasta 4 Bars. Frecuencia hasta 15 Hz.

EL.81-CH: Presión hasta 5 Bars. Frecuencia hasta 20 Hz

**Accesorios incluidos****EL.80-CH:**

- 1 Aplicador SWT
- 1 Transmisor multifocalizado de 9 mm
- 1 Transmisor focalizado de 15 mm
- 1 Transmisor multifocalizado de 15 mm
- 1 Gel 260 ml
- 1 Llave USB 4 Gb

**Accesorios incluidos****EL.81-CH:**

- 1 Aplicador SWT
- 1 Transmisor multifocalizado de 9 mm
- 1 Transmisor focalizado de 15 mm
- 1 Transmisor multifocalizado de 15 mm
- 1 Gel 260 ml
- 1 Llave USB de 4GB
- 1 Kit intercambiable



## EL.81-CH

## EL.80-CH



<b>Línea</b>		
<b>Pantalla</b>	Gráfica Touch Screen a color 8'	Gráfica Touch Screen a color 8'
<b>Presión</b>	1- 5 Bar	1- 4 Bar
<b>Frecuencia</b>	Max 20 Hz	Max 15 Hz
<b>Período de tratamiento</b>	0 - 2 millones de coque	0 - 2 millones de coque
<b>Protocolos memorizados</b>	20	20
<b>Protocolos memorizables en la memoria interna</b>	200	200
<b>Protocolos memorizables en la Llave USB</b>	X	X

## Accesorios en dotación

	EL.81-CH	EL.80-CH
Aplicador SWT Shock Med	X	X
Kit intercambiable	X	-
Transmisor multifocalizado 9 mm	X	X
Transmisor multifocalizado 15 mm	X	X
Transmisor focalizado 15 mm	X	X
Gel 100ml	X	X
Llave USB	X	-
Smart Card	-	X
Manual de uso en DVD	X	X
Cable de alimentación	X	X

## Accesorios opcionales

	CODIGO
Aplicador con trasmisor de 15 mm	ACC1323
Kit intercambiable aplicador	ACC1323/0
Trasmisor focalizado de 15 mm	ACC1273
Trasmisor multifocalizado de 9 mm	ACC1273/0
Trasmisor multifocalizado de 15 mm	ACC1273/1
Llave USB 4GB EME	ACC1317/2/E
Maleta vacía en TNT para transporte	CONT72
Carro con 3 estantes	ACC604

# CARACTERÍSTICAS TÉCNICAS / TECHNICAL SPECIFICATIONS

MODELS	EL.81-CH	EL.80-CH
Main power supply	230 Vca, 50-60 Hz ± 10% 115 Vca, 50-60 Hz ± 10%*	230 Vca, 50-60 Hz ± 10% 115 Vca, 50-60 Hz ± 10%*
Maximum mains power absorption	80 VA	80 VA
Maximum compressor absorption	550 VA	550 VA
Frequency	from 50Hz to 60Hz	from 50Hz to 60Hz
Backlit LCD display	Graphic colour LCD 8" Touch Screen	Graphic colour LCD 6" Touch Screen
Risk class (93/42/CEE)	IIB	IIB
Intensity of electrical discharge	1- 5 Bar max	1- 4 Bar max
Frequency of electrical discharge	1- 20 Hz max	1- 15 Hz max
Treatment period	0 - 2 millions rounds	0 - 2 millions rounds
Stored protocols	20	20
Storable protocols on the internal memory	200	200
Storable protocols on the USB	4 Gb	
Patient's Card	•	•
One-touch mode	•	•
Weight	31 Kg	8 Kg
Size	39 x 91 x 32 cm	39 x 17,5 x 28 cm
*on request		

**SHOCKWAVES UNIT (4 BARS, 15 Hz.) – EL.80-CH****SHOCKWAVES UNIT (5 BARS, 20 Hz.) – EL.81-CH**

20 protocols in memory.

Internal memory with memory capacity of up to 200 programs.

Voltage: 220-240V / 50-60Hz. On request 110-127V / 50-60Hz.

EL.80-CH: Pressure up to max. 4 bars. Frequency up to 15 Hz

EL.81-CH: Pressure up to max. 5 bars. Frequency up to 20 Hz

**Accessories included****EL.80-CH:**

- 1 SWT applicator.
- 1 x 9 mm multi-focused transmitter
- 1 x 15mm focused transmitter
- 1 15 mm multi-focused transmitter
- 1 Gel 260 ml
- 1 4 Gb USB key

**Accessories included****EL.81-CH:**

- 1 SWT applicator.
- 1 x 9 mm multi-focused transmitter
- 1 x 15mm focused transmitter
- 1 15 mm multi-focused transmitter
- 1 Gel 260 ml
- 1 4 Gb USB key
- 1 interchangeable kit.





## Ortotecsa,S.L.

C/Peña de Francia,8-8bis. CP: 28500 Arganda del Rey, Madrid (España)

Tlf: 0034 918 719 132 / Fax: 0034 918 719 323

Email: [info@ortotecsa.com](mailto:info@ortotecsa.com) / Web: [www.ortotecsa.com](http://www.ortotecsa.com)

# Shock waves

	EL.81-CH	EL.80-CH
<b>Line</b>		
<b>Display</b>	Colour Touch Screen 8'	Colour Touch Screen 6'
<b>Pression</b>	Up to 5 Bar	Up to 4 Bar
<b>Frequency</b>	Max 21 Hz	Max 15 Hz
<b>Shots guaranteed</b>	About 2 million	About 2 million
<b>Stored protocols</b>	20	20
<b>Storable protocols on the internal memory</b>	200	200
<b>Storable protocols on the USB</b>	X	X

## Supplied accessories

	EL.81-CH	EL.80-CH
Applicator with 15 mm transmitter	X	X
Interchangeable kit for the applicator	X	-
Multi focused transmitter, 9mm	X	X
Multi focused transmitter, 15mm	X	X
Focused transmitter, 15 mm	X	X
Gel 100ml	X	X
USB key	X	-
Smart Card	-	X
DVD User's Manual	X	X

## Optional accessories

	CODE
Applicator with 15 mm transmitter	ACC1323
Interchangeable kit for the applicator	ACC1323/0
Focused transmitter, 15 mm	ACC1273
Multi focused transmitter, 9mm	ACC1273/0
Multi focused transmitter, 15mm	ACC1273/1
Usb key	ACC1317/2/E
Transport bag in tnt fabric	CONT72
3-shelf trolley	ACC604



# CARACTERÍSTICAS TÉCNICAS / TECHNICAL SPECIFICATIONS

MODELS	EL.81-CH	EL.80-CH
Main power supply	230 Vca, 50-60 Hz ± 10% 115 Vca, 50-60 Hz ± 10%*	230 Vca, 50-60 Hz ± 10% 115 Vca, 50-60 Hz ± 10%*
Maximum mains power absorption	80 VA	80 VA
Maximum compressor absorption	550 VA	550 VA
Frequency	from 50Hz to 60Hz	from 50Hz to 60Hz
Backlit LCD display	Graphic colour LCD 8" Touch Screen	Graphic colour LCD 6" Touch Screen
Risk class (93/42/CEE)	IIB	IIB
Intensity of electrical discharge	1- 5 Bar max	1- 4 Bar max
Frequency of electrical discharge	1- 20 Hz max	1- 15 Hz max
Treatment period	0 - 2 millions rounds	0 - 2 millions rounds
Stored protocols	20	20
Storable protocols on the internal memory	200	200
Storable protocols on the USB	4 Gb	
Patient's Card	•	•
One-touch mode	•	•
Weight	31 Kg	8 Kg
Size	39 x 91 x 32 cm	39 x 17,5 x 28 cm
*on request		